ABUTMENT LINE

IFU_AL0013_ML_REV.C_2015-02-04

PreFace-Abutments

2

DE

PreFace-Abutments

4

EN

Piliers PreFace

6

FR

Monconi PreFace

8

IT

Pilares PreFace

10

ES

Abutments PreFace

12

PT

Hersteller
Manufacturer
Fabricant
Produttore
Fabricante
Fabricante

Medentika® GmbH
Hammweg 8-10
DE-76549 Hügelsheim/Germany
Tel: +49-7229-69912-0
Fax: +49-7229-69912-20
E-Mail: info@medentika.de
www.medentika.de

PreFace-Abutments

B-, C-, D-, E-, F-, H-, I-, K-, L-, N-, R-, S-, T- und Y-Serie

Indikation:

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B-, C-, D-, E-, F-, H-, I-, K-, L-, N-, R-, S-, T and Y-Series

Indication:

For the production of individual titanium/CoCr abutments on implants (using CAD/CAM milling machines). The individual titanium/CoCr abutments in combination with crowns and bridges, for reconstruction of function and aesthetics be made.

The PreFace-Abutments of the B-Series are indicated for Bredent Medical® / SKY® The current implant diameters has to be attended.
The PreFace-Abutments of the C-Series are indicated for Altatec® / Camlog®. The current implant diameters has to be attended.
The PreFace-Abutments of the D-Series are indicated for Altatec® / Conelog®. The current implant diameters has to be attended.
The PreFace-Abutments of the E-Series are indicated for Nobel Biocare® / NobelReplace® Tapered. The current implant diameters has to be attended.
The PreFace-Abutments of the F-Series are indicated for Nobel Biocare® / NobelActive® - NobelReplace® Conical. The current implant diameters has to be attended.
The PreFace-Abutments of the H-Series are indicated for BIOMET 3i® / Certain®. The current implant diameters has to be attended.
The PreFace-Abutments of the I-Series are indicated for BIOMET 3i® / External Hex. The current implant diameters has to be attended.
The PreFace-Abutments of the K-Series are indicated for Nobel Biocare® / Brånemark® The current implant diameters has to be attended.
The PreFace-Abutments of the L-Series are indicated for Straumann® / Bone Level®. The current implant diameters has to be attended.
The PreFace-Abutments of the N-Series are indicated for Straumann® / Soft Tissue Level® The current implant diameters has to be attended.
The PreFace-Abutments of the R-Series are indicated for Zimmer Dental® / Tapered Screw-Vent®, MIS® / SEVEN® Internal Hex, BioHorizons® / Tapered Internal - Tapered Internal Plus - Tapered Tissue Level. The current implant diameters has to be attended.
The PreFace-Abutments of the S-Series are indicated for DENTSPLY Implants® / ASTRA TECH OsseoSpeed® TX The current implant diameters has to be attended.
The PreFace-Abutments of the T-Series are indicated for DENTSPLY Implants® / XiVE® S The current implant diameters has to be attended.
The PreFace-Abutments of the Y-Series are indicated for DENTSPLY Implants® / ANKYLOS® The current implant diameters has to be attended.

The PreFace-Abutments are provided with screw included.

Composition:

PreFace-Abutment titanium:

Ti6Al4V, medical grade 5, ASTM 136

PreFace-Abutment CoCr:

CoCr28Mo, ASTM F899, WAK 14,1

Abutment screw: Ti6Al4V, medical grade 5, ASTM 136

Contraindication:

The PreFace-Abutments of each series can only be combined with the corresponding compatible implant system, e.g. the E-series only in combination with Nobel Biocare® / NobelReplace® Tapered. No inappropriate in connection geometry abutments should be used. Any post-processing at the connection geometry to the implant may result in fitting inaccuracies prohibiting further use. All PreFace-Abutments are intended for single use. There should be no other than that of the manufacturer's recommended torque be used.

Ncm	Abutment
15	F-Series 3,0; S-Series 3,0
20	C-Series D-Series H-Series
25	B-Series; S-Series 3,5/4,0 4,5/5,0; T-Series; Y-Series
30	R-Series
35	I-Series; K-Series; L-Series; N-Series; E-Series; F-Series NP/RP

C€0483

In general:

PreFace abutments are processed in a suitable milling system using original PreFace abutment holders. The PreFace abutment holder is for gripping the PreFace abutment in the correct position in the milling machine in order to process it as required.

PreFace abutments are gripped face-side first in their PreFace abutment holders ensuring that they are fully seated. PreFace abutment holders are available for many milling machines (see current PreFace catalogue).

Designing:

The PreFace-Abutment may at its outer form be individually adapted to the anatomical condition. The production of the individual PreFace-Abutment done digitally with a suitable CAD software. Basic rules should be observed. The design of the emergence profile should be done gently. A greater widening of the gingiva, was formed as by healing caps should always be done after consultation with the doctor treated. The gingiva can be expanded gradually by custom-made temporary abutments. In general the following applies: Place the circumferential shoulder in the labial region slightly subgingivally and follow the gingival contour in the oral region. This allows the cement between the copings and abutment to be easily removed and the junction is not visible from the labial aspect. The abutment thickness should not be less than 0.5 mm (screw hole to outer abutment surface).

The diameter of the emergence profile should 13 mm, the height of the emergence profile should 6 mm not exceed. The total length of the PreFace-Abutment should 15 mm, the maximum angulation of the PreFace-Abutment should 30° not exceed. The taper of the PreFace-Abutment should not exceed 6°. The PreFace-Abutment should be rounded occlusal and not leak pointed. Sharp edges should be avoided.

Handling:

If further steps are required following processing of the PreFace abutment in the CAD/CAM milling machine, please heed the following:

Appropriate milling machines, in a faultless condition, must be used under slight pressure to machine the PreFace-Abutments. The minimum wall thickness of 0.5 mm may not be undershot. Ridges and edges must be avoided.

We recommend that you affix the PreFace-Abutments in laboratory implants to safeguard the connection geometry.

Veneering:

If an PreFace-Abutment is to be veneered directly, the veneering is produced using suitable veneering materials in accordance with the relevant manufacturer's instructions.

Scanbody:

For detection of the precise implant position during scanning, use the scanbody 2. Generation. This must be selected to be compatible with the original implant system and diameter according to the titanium base range to be used.

Safety notice:

Metal dust is harmful to health. When machining and sandblasting use a suction with a fine dust filter that is usual in practice and wear protective goggles as well as a face mask.

Side effects:

Allergies or sensibilities in connection with the alloy cannot be ruled out in very rare instances.

Interactions:

Various alloy types in the same mouth cavity can lead to galvanic reactions in the event of occlusal or proximal contact.

Notice:

The information conveyed orally, in writing or in practical seminars is based on tests and experience and can therefore only be considered to be standard values. Our products are subject to constant further development. In this connection we reserve the right to make product changes in relation to the design and composition.

Cleaning and sterilization:

The PreFace-Abutments of the B, C, D, E, F, H, I, K, L, N, R, S, T, Y series are delivered in a non-sterile state in the appropriate packaging. The products has to be cleaned and sterilized according to following details:

Cleaning:

Automatic cleaning procedure based on Vario TD program:

- Cold water rinse 1 min.
- Cleaning at 55°C/131°F (± 2°C/35,6°F) for 5 min. with cleaning agent neodisher® MediClean (0.5% v/v)
- Neutralization with cold water for 2 min.
- Cold water rinse 1 min.

PreFace-Abutments

B-, C-, D-, E-, F-, H-, I-, K-, L-, N-, R-, S-, T and Y-Series

C €0483

Sterilization:

Seal in a suitable autoclave foil.

Steam sterilization with 3x fractionated pre-vacuum:

132°C (270°F) for 4 minutes, drying time: 20 minutes

Warning:

The product has to be inspected prior usage. The packing must be sealed without any visible damage. The following descriptions are not sufficient for the immediate use of the product. Dental skills and prior instruction of how to use the product are at any rate required.

Signs and Symbols:

For single use only

LOT number

Reference number

Manufacturer

Attend instruction for use

Non sterile

Products indicated with ® are registered brand names of the manufacturers.

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